

Clarisonic Opal is hydrating and gentle

Study Objective:

Assess gentleness of the Clarisonic Opal Sonic Infusion System for twice daily use

Methodology:

- 30 subjects participated in this home use study
- A board certified dermatologist and board certified ophthalmologist examined the skin and eyes at baseline and after two weeks of use
- Participants cleansed their face using their normal skin care routine
- After cleansing, participants used Clarisonic Opal with the Anti-Aging Sea Serum for 30 seconds around each eye area twice a day
- Non-invasive skin barrier measurements (Transepidermal Water Loss, Erythema, Hydration) on the crow's feet area of both eyes were also taken at baseline and after two weeks of use

Results:

Improvement in the skin barrier (increased skin hydration: $p < 0.01$ and reduced Transepidermal Water Loss: $p < 0.01$) was reported after 2 weeks use of the Clarisonic Opal and Anti-Aging Sea Serum

Dermatologist and ophthalmologist assessments confirmed that Clarisonic Opal is gentle for daily use

